



## Application Guide

### Abri-San (Shaped) – Upright position



1. Fold the pad lengthways with the backsheet turned outwards. Do not touch the inside of the pad.
2. Ensuring the user is properly supported pass the pad from front to back. Position the upper rear end of the pad at the coccyx.

3. Fan out the back of the pad horizontally to create a channel shape between the legs.



4. Fan out pad horizontally at the front.



5. Pull the fitting pants into position. If possible, the pant legs should be placed on the backsheet edges of the pad.



6. Adjust the fitting pants at the front and make sure the edges of the pad ease into the groin area.